

**PREG-TONE II<sup>®</sup>**  
*plus*  
“Ultrasound for Ultra Profits”<sup>®</sup>

**fast**

Testing can be completed in  
less than 30 seconds.

**safe**

All testing done externally —  
no harm to sow or litter.

**accurate**

Accuracy proven by many  
thousands of producers.

**DISTRIBUTED BY:**

Renco Corporation  
116 Third Ave. North  
Minneapolis, Minnesota 55401, USA

**Tel:**

612-338-6124

**Toll Free:**

1-800-359-8181

**Fax:**

612-333-9026

**Website:**

[www.renecorp.com](http://www.renecorp.com)

**Email:**

[sales@renecorp.com](mailto:sales@renecorp.com)



FORM 1000

“Ultrasound for Ultra Profits”<sup>®</sup>

©2000 RENCO CORPORATION. ALL RIGHTS RESERVED.



**PREG-TONE II<sup>®</sup>**  
*plus*

**ANIMAL PREGNANCY TESTER**



**DETECTS PREGNANCY IN SWINE**

✓ *fast*   ✓ *accurate*   ✓ *efficient*

# PREG-TONE II<sup>®</sup> plus

AVAILABLE WITH OPTIONAL EARPHONE

*Easy one-hand operation*



*The Preg-Tone II Plus is very simple to use.*

**Step one:** Apply contact fluid, such as mineral or cooking oil, to animal in area shown.

**Step two:** Turn the instrument on by pressing the patented case with your thumb.

**Step three:** Press instrument to the skin and scan interior of the animal.

## DURABLE

- *High-impact plastic*
- *Acid resistant*
- *Will last for years*

## INEXPENSIVE

- *9 volt replaceable battery*
- *Battery will last for months*

## EASY TO USE

- *One hand operation*
- *Use any oil for skin contact*
- *No lights to watch*

## INDICATORS

- *Loud Sounder*
- *Intermittent tone indicates good skin contact*
- *Continuous tone indicates pregnancy*

## OPTIONAL

- *Earphone attachment available for extra noisy conditions*

## QUALITY

- *Renco has over 30 years experience testing and manufacturing ultrasound equipment.*

